

The Problem: Current Standards of Urinary Incontinence Care

The wrong product for urinary incontinence can significantly influence a person's body image and self-esteem, decrease quality of life, and lead to depression and isolation.

Did you know that over 18%¹ or almost twice as many men with military service experience urinary incontinence?

Condom Catheters



- Condom Catheters can expose skin to urine and moisture
- Multiple sizes/types - all can cause **skin inflammation & infections**
- Wearing one too tightly can cause reduced blood flow and tissue damage
- Conversely, wearing one too loosely can cause “pop-offs” and leaks.
- Mechanical factors such as friction between the skin and any type of material contributes to **skin damage**².

Diapers/Pads



- Exposure to high moisture levels for prolonged periods can be detrimental to the skin, leading to inflammation and erosion of the skin.²
- Frequently causes skin inflammation sores and ulcers
- Urine smells as it accumulates
- Can undermine personal dignity
- Expensive, not covered by insurance, and can cost patients up to \$300 a month.

The Solution

Men's Liberty is the first truly innovative device for men's urinary incontinence that keeps men dry 24/7. It features a completely external design, a skin friendly hydrocolloid adhesive and most men get a secure seal for up to 24 hours, with more than 5 million have been sold to date. Men's Liberty is the discreet, dignified and dependable option for managing urinary incontinence.

Patients Men's Liberty Can Help

- Geriatric Incontinence
- Stress Related Incontinence
- Obesity
- Uncircumcised Anatomy
- Retracted Anatomy
- Limited Dexterity
- Diabetes Mellitus
- Small Anatomy
- Neuromuscular syndromes
- Orthopaedic surgery of the pelvis with a normal voiding pattern
- Post surgical incontinence
- Limited Mobility/ Paralyzed and many more

An External Catheter that works

The Men's Liberty™ External Catheter is a one size fits most, premier alternative to condom catheters, diapers & pads, a noninvasive solution to male urinary incontinence delivering superior performance, leak protection, and patient comfort



One size fits most

Fits nearly all men and is suitable for circumcised and uncircumcised patients



Safe on skin

Leakproof design keeps patients dry while proprietary hydrocolloid adhesive provides a strong, gentle adhesion



Extended wear time

Up to 24 hour wear time (after first application of 6–12 hours)³. The hydrocolloid turns white when it needs to be removed



Available through the VA CMOP program.

CMOP # XT660 Item # BOM23406H

Men's Liberty is manufactured by BioDerm, Inc. in Largo, FL and has been approved for billing under HCPCS Code A4326. Men's Liberty works with a network of distribution partners in the US and internationally to bring Men's Liberty straight to men's door steps.

7-petal faceplate made of gentle hydrocolloid adheres directly to patient's skin preventing leakage

Hydrocolloid faceplate seal provides a second layer of leak protection

Gravitational vent allows air, ushering urine away from skin to prevent pooling and skin breakdown

Urine Collection Pouch is discrete and undetectable under the clothing.

Drain Plug allows voiding without removing the product. Can connect to standard bed or leg bag for more capacity.

References: ¹ Vaughan CP, Johnson TM 2nd, Goode PS, Redden DT, Burgio KL, Markland AD. Military exposure and urinary incontinence among American men. J Urol. 2014 Jan; 191(1):125-9. doi: 10.1016/j.juro.2013.07.016. Epub 2013 Jul 16. PMID: 23871759. ² Surber C et al., 2015; ³ Beekman D et al., 2016.

For Questions & Free Sample Requests:
www.MensLiberty.com/Veterans